Uputa za Cystatin C imunotest na sustavima Beckman Coulter® AU (AU5800, AU680, AU480, AU2700)

REF B08179

Namjena

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Instrument Settings for Cystatin C Immunoassay

Primjene na serumu i plazmi uređaja Cystatin C AU2700

Sistemske reagens: B08179 ID reagens: 228

Specific Test Parameters									
General		LIH	ISE	Range					
Test Name:	CysC	<	>	Type:	Serum	Operation:	Yes		
Sample:	Volume	2	μL	Dilution	0	μL	Pre-Dilution Rate:	1	
Reagents:	R1 Volume	150	μL	Dilution	0	μL	Min OD		Max OD
	R2 Volume	30	μL	Dilution	10	μL	L		H
Wavelength:	Pri.	540	nm	Sec.		nm	Reagent OD limit:		
Method:	End Point			First L	-2.0	First H	2.0	Last L	Last H
Reaction slope:		+		Dynamic Range:	L	*	H	*	
Measuring Point 1:	First	13		Last	27		Correlation Factor:	A	1.00
Measuring Point 2:	First			Last			B	0.00	
Linearity:			%	On-board stability period:	60				
No Lag Time:									

Specific Test Parameters									
General		LIH	ISE	Range					
Test Name:	CysC	<	>	Type:	Serum				
Value/Flag:	#	Level L:	#	Level H:	#				
Normal Ranges:	Age L		Age H		L		H		
	Sex	Year	Month	Year	Month				
o 1.	#	#	#	#	#	#	#	#	#
o 2.	#	#	#	#	#	#	#	#	#
o 3.	#	#	#	#	#	#	#	#	#
o 4.	#	#	#	#	#	#	#	#	#
o 5.	#	#	#	#	#	#	#	#	#
o 6.	#	#	#	#	#	#	#	#	#
7.	None Selected								
8.	Out of Range		L	H					
Panic Value:	#	#	#	Unit:	mg/L	Decimal places:	#		

Calibration Specific									
General		ISE							
Test Name:	CysC	<	>	Type:	Serum				
Calibration Type:	6AB	Formula:	Spline	Counts:	#	Process:	CONC		
Point 1:	Cal. No.	OD	CONC	Factor/OD-L	Factor/OD-H				
Point 2:	1		**						
Point 3:	2		**						
Point 4:	3		**						
Point 5:	4		**						
Point 6:	5		**						
Point 7:	6		**						
1-Point Cal. Point:		o	with CONC-0	Slope Check:	+	Advanced Calibration:	#		
MB Type Factor:		Calibration Stability Period:	28						

Definira korisnik

** Specifično za seriju, vidjeti list s analitičkim vrijednostima uključen u komplet s kalibratorom

Primjene na serumu i plazmi uređaja Cystatin C AU680/AU480:

Sistemske reagensi: B08179 ID reagensi: 228

Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation: <input type="text" value="Yes"/>					
Sample Volume	<input type="text" value="2"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD	Max.OD
Rgt. Volume	R1(R1-1) <input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit	
				First	Low
				Last	Low
					High
					High
R2(R2-1)	<input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Dynamic Range Low	
					High
Common Rgt. Type	<input type="text"/>	Name	<input type="text"/>	Correlation Factor A	<input type="text" value="1.00"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec.	<input type="text"/>	Factor for Maker A	<input type="text" value="B"/>
Method	<input type="text" value="End Point"/>				
Reaction Slope	<input type="text" value="+"/> ∇	Onboard Stability Period <input type="text" value="60"/> Day <input type="text"/> Hour			
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check <input type="text"/>	
Measuring Point2 First	<input type="text"/>	Lipemia <input type="text"/>			
Linearity Limit	<input type="text"/>	Icterus <input type="text"/>			
Lag Time Check	<input type="text"/>	Hemolysis <input type="text"/>			

Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>					
Value/Flag:	<input type="text" value="#"/>	Low	<input type="text"/>	High	<input type="text"/>
Level	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Specific Ranges:	From	To	Panic Value		
	Sex	Year	Month	Year	Month
o 1.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
o 2.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
o 3.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
o 4.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
o 5.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
o 6.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
o 7.	No demographics				<input type="text"/>
o 8.	Not within expected values				<input type="text"/>
Unit	<input type="text" value="mg/L"/>	Decimal Places	<input type="text" value="#"/>		

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> <input type="checkbox"/> Use Serum Cal.					
Calibration Type:	<input type="text" value="6AB"/>	Formula:	<input type="text" value="Spline"/>	Counts:	<input type="text" value="#"/>
<Calibrator Parameters>					
Calibrator	OD	Conc	Low	High	Slope Check
Point 1:	<input type="text" value="1"/>	<input type="text" value="**"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="+"/> ∇
Point 2:	<input type="text" value="2"/>	<input type="text" value="**"/>	<input type="text"/>	<input type="text"/>	Allowance Range Check
Point 3:	<input type="text" value="3"/>	<input type="text" value="**"/>	<input type="text"/>	<input type="text"/>	<input type="checkbox"/> Reagent Blank
Point 4:	<input type="text" value="4"/>	<input type="text" value="**"/>	<input type="text"/>	<input type="text"/>	<input type="checkbox"/> Calibration
Point 5:	<input type="text" value="5"/>	<input type="text" value="**"/>	<input type="text"/>	<input type="text"/>	Advanced Calibration
Point 6:	<input type="text" value="6"/>	<input type="text" value="**"/>	<input type="text"/>	<input type="text"/>	Operation
Point 7:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="#"/> ∇
Point 8:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Interval (RB/ACAL)
Point 9:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="#"/> ∇
Point 10:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
<Point Cal. For	No. of Correction Points	<input type="text"/>	Use Master Curve	<input type="checkbox"/>	<input type="checkbox"/> Lot Calibration
Master Curve>					
Calibrator	OD	Conc	Low	High	Stability
Point-1	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Reagent Blank
Point-2	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Calibration
MB Type Factor:	<input type="text"/>	1-Point Calibration Point	<input type="checkbox"/>	<input type="checkbox"/>	<input type="text" value="28"/> Day <input type="text"/> Hour
					<input type="text" value="28"/> Day <input type="text"/> Hour
					<input type="checkbox"/> with Conc-0

Definira korisnik

** Specifično za seriju, vidjeti list s analitičkim vrijednostima uključen u komplet s kalibratorom

Primjene na serumu i plazmi uređaja Cystatin C AU5800

Sistemske reagens: B08179 ID reagens: 228

Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation <input type="text" value="Yes"/>					
Sample Volume	<input type="text" value="2"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>	Diluent Bottle	<input type="text" value="Outside"/>	Min.OD	Max.OD
Rgt. Volume	<input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit	
R1-2	<input type="text" value=""/>	Dilution	<input type="text" value=""/>	First	Low
				Last	Low
R2(R2-1)	<input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Dynamic Range Low	High
Common Rgt. Type	<input type="text" value=""/>	Name	<input type="text" value=""/>	Correlation Factor A	B
Wavelength	Pri <input type="text" value="540"/> nm	Sec.	<input type="text" value=""/>	Factor for Maker A	B
Method	<input type="text" value="End Point"/>				
Reaction Slope	<input type="text" value="+"/>				
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	Onboard Stability Period	<input type="text" value="60"/> Day <input type="text" value=""/> Hour
Measuring Point2 First	<input type="text" value=""/>	Last	<input type="text" value=""/>	LIH Influence Check	<input type="text" value=""/>
Linearity Limit	<input type="text" value=""/>				
Lag Time Check	<input type="text" value=""/>				
				Lipemia	<input type="text" value=""/>
				Icterus	<input type="text" value=""/>
				Hemolysis	<input type="text" value=""/>

Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>					
Value/Flag: <input type="text" value="#"/>					
Level					
Specific Ranges:		From	To	Low	High
	Sex	Year	Month	Year	Month
o 1.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
o 2.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
o 3.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
o 4.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
o 5.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
o 6.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
7.	Standard demographics				
8.	Not within expected values				
Panic Value	Low	<input type="text" value=""/>	High	<input type="text" value=""/>	Unit <input type="text" value="mg/L"/> Decimal Places <input type="text" value="#"/>

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Cuvette: <input type="text" value=""/>					
<input type="checkbox"/> Use Serum Cal.					
Calibration Type: <input type="text" value="6AB"/> Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="#"/>					
<Calibrator Parameters>					
Point 1:	Calibrator	OD	Conc	Factor Range	Slope Check
Point 2:	<input type="text" value="1"/>	<input type="text" value=""/>	<input type="text" value="**"/>	Low	High
Point 3:	<input type="text" value="2"/>	<input type="text" value=""/>	<input type="text" value="**"/>		
Point 4:	<input type="text" value="3"/>	<input type="text" value=""/>	<input type="text" value="**"/>		
Point 5:	<input type="text" value="4"/>	<input type="text" value=""/>	<input type="text" value="**"/>		
Point 6:	<input type="text" value="5"/>	<input type="text" value=""/>	<input type="text" value="**"/>		
Point 7:	<input type="text" value="6"/>	<input type="text" value=""/>	<input type="text" value="**"/>		
Point 8:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>		
Point 9:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>		
Point 10:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>		
Allowance Range Check					
<input type="checkbox"/> Reagent Blank <input type="text" value=""/>					
<input type="checkbox"/> Calibration <input type="text" value=""/>					
Advanced Calibration					
Operation <input type="text" value="#"/>					
Interval (RB/ACAL) <input type="text" value="#"/>					
<Point Cal. For No. of Correction Points <input type="text" value=""/> Use Master Curve <input type="checkbox"/> Lot Calibration <input type="checkbox"/>					
Master Curve>					
Point-1	Calibrator	OD	Conc	Low	High
Point-2	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>		
Stability					
Reagent Blank <input type="text" value="28"/> Day <input type="text" value="0"/> Hour					
Calibration <input type="text" value="28"/> Day <input type="text" value="0"/> Hour					
MB Type Factor: <input type="text" value=""/> 1-Point Calibration Point <input type="checkbox"/> with Conc-0 <input type="checkbox"/>					

Definira korisnik

** Specifično za seriju, vidjeti list s analitičkim vrijednostima uključen u komplet s kalibratorom