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Calibration Specific											
General		ISE									
Test Name:		CysC		< >		Type:		Serum			
Calibration Type:		6AB		Formula:		Spline		Counts: #		Process: CONC	
Point 1:		Cal. No. 1		OD		CONC **		Factor/OD-L		Factor/OD-H	
Point 2:		2				**					
Point 3:		3				**					
Point 4:		4				**					
Point 5:		5				**					
Point 6:		6				**					
Point 7:						**					
1-Point Cal. Point:		<input type="checkbox"/>		with CONC-0		Slope Check: +		Advanced Calibration: #			
MB Type Factor:				Calibration Stability Period:		28					

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Aplikace cystatin C AU680/AU480 pro sérum a plasmu

Systémová reagenzie: B08179 ID reagenzie: 228

Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < >		Type: <input type="text" value="Serum"/>		Operation: <input type="text" value="Yes"/>	
Sample Volume	<input type="text" value="2"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD	<input type="text"/>
Rgt. Volume R1(R1-1)	<input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit	<input type="text"/>
				First	Low: <input type="text" value="-2.0"/> High: <input type="text" value="2.0"/>
				Last	Low: <input type="text"/> High: <input type="text"/>
R2(R2-1)	<input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Dynamic Range Low	0.4 High: <input type="text" value="7.8"/>
Common Rgt. Type	<input type="text"/>	Name	<input type="text"/>	Correlation Factor A	1.00 B: <input type="text" value="0.00"/>
Wavelength	Pri: <input type="text" value="540"/> nm	Sec.	<input type="text"/>	Factor for Maker A	<input type="text"/>
Method	<input type="text" value="End Point"/>				
Reaction Slope	<input type="text" value="+"/> ∇			Onboard Stability Period	<input type="text" value="60"/> Day <input type="text"/> Hour
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check	<input type="text"/>
Measuring Point2 First	<input type="text"/>	Last	<input type="text"/>	Lipemia	<input type="text"/>
Linearity Limit	<input type="text"/>			Icterus	<input type="text"/>
Lag Time Check	<input type="text"/>			Hemolysis	<input type="text"/>

Parameters		Specific Test Parameters						
General	LIH	ISE	Calculated Test	Range				
Test Name: <input type="text" value="CysC"/> < >		Type: <input type="text" value="Serum"/>						
Value/Flag:	<input type="text" value="#"/>	Low	<input type="text"/>	High	<input type="text"/>			
Level	<input type="text" value="#"/>							
Specific Ranges:	From	To			Panic Value			
	Sex	Year	Month	Year	Month	Low	High	
<input type="checkbox"/> 1.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
<input type="checkbox"/> 2.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
<input type="checkbox"/> 3.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
<input type="checkbox"/> 4.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
<input type="checkbox"/> 5.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
<input type="checkbox"/> 6.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
7.	No demographics						<input type="text" value="#"/>	<input type="text" value="#"/>
8.	Not within expected values						<input type="text" value="#"/>	<input type="text" value="#"/>
Unit	<input type="text" value="mg/L"/>		Decimal Places	<input type="text" value="#"/>				

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="CysC"/> < >		Type: <input type="text" value="Serum"/>		<input type="checkbox"/> Use Serum Cal.	
Calibration Type: <input type="text" value="6AB"/>		Formula: <input type="text" value="Spline"/>		Counts: <input type="text" value="#"/>	
<Calibrator Parameters>					
Calibrator	OD	Conc	Factor Range	Low	High
Point 1:	<input type="text" value="1"/>	<input type="text" value="**"/>			
Point 2:	<input type="text" value="2"/>	<input type="text" value="**"/>			
Point 3:	<input type="text" value="3"/>	<input type="text" value="**"/>			
Point 4:	<input type="text" value="4"/>	<input type="text" value="**"/>			
Point 5:	<input type="text" value="5"/>	<input type="text" value="**"/>			
Point 6:	<input type="text" value="6"/>	<input type="text" value="**"/>			
Point 7:	<input type="text"/>	<input type="text"/>			
Point 8:	<input type="text"/>	<input type="text"/>			
Point 9:	<input type="text"/>	<input type="text"/>			
Point 10:	<input type="text"/>	<input type="text"/>			
Slope Check: <input type="text" value="+"/>					
Allowance Range Check					
<input type="checkbox"/> Reagent Blank <input type="text"/>					
<input type="checkbox"/> Calibration <input type="text"/>					
Advanced Calibration Operation: <input type="text" value="#"/>					
Interval (RB/ACAL): <input type="text" value="#"/>					
<Point Cal. For Master Curve>					
Calibrator	OD	Conc	OD Range	Low	High
Point-1	<input type="text"/>	<input type="text"/>			
Point-2	<input type="text"/>	<input type="text"/>			
Stability: Reagent Blank <input type="text" value="28"/> Day <input type="text"/> Hour					
Calibration: <input type="text" value="28"/> Day <input type="text"/> Hour					
MB Type Factor: <input type="text"/> 1-Point Calibration Point <input type="text"/>					
<input type="checkbox"/> with Conc-0					

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ID reagenčie: 228

General	LIH	ISE	Calculated Test	Range
Test Name: CYS1G Test No Type: Serum Operation: Yes				
Sample Volume: <input type="text" value="2.0"/> µL	Dilution: <input type="text" value="0"/> µL	OD Limit	Min. OD: <input type="text"/>	Max OD: <input type="text"/>
Pre-Dilution Rate: <input type="text" value="1"/>		Reagent OD Limit	1 st Low: <input type="text" value="-2.0000"/> High: <input type="text" value="2.0000"/>	
Reagent Volume R1 (R1-1): <input type="text" value="150"/> µL	Dilution: <input type="text" value="0"/> µL	Last Low: <input type="text" value="-2.0000"/> High: <input type="text" value="2.0000"/>	Analytical Measuring Range Low: <input type="text" value="0.4"/> High: <input type="text" value="7.8"/>	
R1-2: <input type="text"/> µL	Dilution: <input type="text"/> µL	Correlation Factor A: <input type="text" value="1"/> B: <input type="text" value="0"/>		
R2 (R2-1): <input type="text" value="30"/> µL	Dilution: <input type="text" value="10"/> µL	Manufacturer Factor A: <input type="text" value="1"/> B: <input type="text" value="0"/>		
Common Reagent Type: <input type="text" value="None"/>	Name: <input type="text" value="None"/>	Onboard Stability Period: <input type="text" value="60"/> Day <input type="text" value="0"/> Hour	LIH Influence Check: <input type="text" value="No"/>	
Wavelength Pri: <input type="text" value="540"/> nm	Sec: <input type="text" value="None"/> nm	Lipemia: <input type="text"/>	Icterus: <input type="text"/>	
Method: <input type="text" value="END"/>		Hemolysis: <input type="text"/>		
Reaction Slope: <input type="text" value="+"/> ▼				
Measuring Point-1 1st: <input type="text" value="13"/>	Last: <input type="text" value="27"/>			
Measuring Point-2 1st: <input type="text"/>	Last: <input type="text"/>			
Linearity Limit: <input type="text"/> %				
Lag Time Check: <input type="text"/>				

General	LIH	ISE	Calculated Test	Range
Test Name: CYS1G Test No Type: Serum				
Value/Flag: <input type="text" value="Value"/>	Level: Low	<input type="text" value="-99999.99"/>	High	<input type="text" value="99999.99"/>
Specific Ranges				
<input type="checkbox"/> 1: Sex #	From Year #	Month #	To Year #	Month #
<input type="checkbox"/> 2: #	#	#	#	#
<input type="checkbox"/> 3: #	#	#	#	#
<input type="checkbox"/> 4: #	#	#	#	#
<input type="checkbox"/> 5: #	#	#	#	#
<input type="checkbox"/> 6: #	#	#	#	#
7: Standard demographics				
8: Not within expected values				
Critical Limits Low: <input type="text" value="#"/>	High: <input type="text" value="#"/>	Unit: <input type="text" value="mg/L"/>	Select: <input type="text"/>	Decimal Places: <input type="text" value="2"/>

Calibrators	General	ISE			
Test Name: CYS1G Type: Serum					
<input type="checkbox"/> Use Serum Cal.					
Calibration Type: <input type="text" value="6AB"/>	Formula: <input type="text" value="Spline"/>	Counts: <input type="text" value="2"/>	Slope Check: <input type="text" value="+"/> ▼		
<Calibrator Parameters>					
Point	Calibrator	OD	Conc	Range	
				Low	High
Point-1	CYSC Calibrator Level 1		†	-2.0000	2.0000
Point-2	CYSC Calibrator Level 2		†	-2.0000	2.0000
Point-3	CYSC Calibrator Level 3		†	-2.0000	2.0000
Point-4	CYSC Calibrator Level 4		†	-2.0000	2.0000
Point-5	CYSC Calibrator Level 5		†	-2.0000	2.0000
Point-6	CYSC Calibrator Level 6		†	-2.0000	2.0000
Point-7					
MB Type Factor: <input type="text"/>	1-Point Calibration Point: <input type="text" value="None"/>	<input type="checkbox"/> with Conc-0		Interval (RB): <input type="text"/>	
				Interval (ACAL): <input type="text"/>	
		Stability		Reagent Blank: <input type="text" value="28"/> Day <input type="text" value="0"/> Hour	
				Calibration: <input type="text" value="28"/> Day <input type="text" value="0"/> Hour	

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