# Cystatine C

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## Cystatine C AU2700/AU5400 Serum en Plasma Toepassing

Systemreagens: B08179 Reagens ID: 228

Specific Test Parameters	
General	LIH ISE Range
Test Name:	CysC < > Type: Serum Operation: Yes
Sample:	Volume: 2 µL Dilution: 0 µL Pre-Dilution Rate: 1
Reagents:	R1 Volume: 150 µL Dilution: 0 µL Min OD Max OD
	R2 Volume: 30 µL Dilution: 10 µL L H
Wavelength:	Pri: 540 Sec: Reagent OD limit: First L -2.0 First H 2.0
Method:	End Point Last L Last H
Reaction slope:	+ Dynamic Range: L 0.4 H 7.8
Measuring Point 1:	First 13 Last 27 Correlation Factor: A 1.00 B 0.00
Measuring Point 2:	First Last
Linearity:	% On-board stability period: 60
No Lag Time:	

Specific Test Parameters	
General	LIH ISE Range
Test Name:	CysC < > Type: Serum
Value/Flag:	# Level L: # Level H: #
Normal Ranges:	Age L Age H L H
	Sex Year Month Year Month L H
<input type="checkbox"/>	1. # # # # # # # #
<input type="checkbox"/>	2. # # # # # # # #
<input type="checkbox"/>	3. # # # # # # # #
<input type="checkbox"/>	4. # # # # # # # #
<input type="checkbox"/>	5. # # # # # # # #
<input type="checkbox"/>	6. # # # # # # # #
	7. None Selected
	8. Out of Range L H
Panic Value:	# # Unit: mg/L Decimal places: #

Calibration Specific	
General	ISE
Test Name:	CysC < > Type: Serum
Calibration Type:	6AB Formula: Spline Counts: # Process: CONC
Point 1:	Cal. No. 1 OD CONC ** Factor/OD-L Factor/OD-H
Point 2:	2 OD CONC ** Factor/OD-L Factor/OD-H
Point 3:	3 OD CONC ** Factor/OD-L Factor/OD-H
Point 4:	4 OD CONC ** Factor/OD-L Factor/OD-H
Point 5:	5 OD CONC ** Factor/OD-L Factor/OD-H
Point 6:	6 OD CONC ** Factor/OD-L Factor/OD-H
Point 7:	7 OD CONC ** Factor/OD-L Factor/OD-H
1-Point Cal. Point:	<input type="checkbox"/> with CONC=0 Slope Check: + Advanced Calibration: #
MB Type Factor:	Calibration Stability Period: 28

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Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> <input type="button" value="&lt;"/> <input type="button" value="&gt;"/>		Type: <input type="text" value="Serum"/>			
Value/Flag: <input type="text" value="#"/>					
Level		Level	Low: <input type="text" value="#"/>	High: <input type="text" value="#"/>	
Specific Ranges:		From	To	Low	High
<input type="checkbox"/> 1.	Sex: <input type="text" value="#"/>	Year: <input type="text" value="#"/>	Month: <input type="text" value="#"/>	Year: <input type="text" value="#"/>	Month: <input type="text" value="#"/>
<input type="checkbox"/> 2.	Year: <input type="text" value="#"/>	Month: <input type="text" value="#"/>	Year: <input type="text" value="#"/>	Month: <input type="text" value="#"/>	Year: <input type="text" value="#"/>
<input type="checkbox"/> 3.	Month: <input type="text" value="#"/>	Year: <input type="text" value="#"/>	Month: <input type="text" value="#"/>	Year: <input type="text" value="#"/>	Month: <input type="text" value="#"/>
<input type="checkbox"/> 4.	Year: <input type="text" value="#"/>	Month: <input type="text" value="#"/>	Year: <input type="text" value="#"/>	Month: <input type="text" value="#"/>	Year: <input type="text" value="#"/>
<input type="checkbox"/> 5.	Month: <input type="text" value="#"/>	Year: <input type="text" value="#"/>	Month: <input type="text" value="#"/>	Year: <input type="text" value="#"/>	Month: <input type="text" value="#"/>
<input type="checkbox"/> 6.	Year: <input type="text" value="#"/>	Month: <input type="text" value="#"/>	Year: <input type="text" value="#"/>	Month: <input type="text" value="#"/>	Year: <input type="text" value="#"/>
7. Standard demographics					
8. Not within expected values					
Panic Value	Low: <input type="text"/>	High: <input type="text"/>	Unit: <input type="text" value="mg/L"/>	Decimal Places: <input type="text" value="#"/>	

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="CysC"/> <input type="button" value="&lt;"/> <input type="button" value="&gt;"/>		Type: <input type="text" value="Serum"/>	Cuvette: <input type="text"/>		
<input type="checkbox"/> Use Serum Cal.					
Calibration Type: <input type="text" value="6AB"/>		Formula: <input type="text" value="Spline"/>	Counts: <input type="text" value="#"/>		
<Calibrator Parameters>					
Point 1:	Calibrator: <input type="text" value="1"/>	OD: <input type="text"/>	Conc: <input type="text" value="**"/>	Factor Range Low: <input type="text"/>	Factor Range High: <input type="text"/>
Point 2:	Calibrator: <input type="text" value="2"/>	OD: <input type="text"/>	Conc: <input type="text" value="**"/>		
Point 3:	Calibrator: <input type="text" value="3"/>	OD: <input type="text"/>	Conc: <input type="text" value="**"/>		
Point 4:	Calibrator: <input type="text" value="4"/>	OD: <input type="text"/>	Conc: <input type="text" value="**"/>		
Point 5:	Calibrator: <input type="text" value="5"/>	OD: <input type="text"/>	Conc: <input type="text" value="**"/>		
Point 6:	Calibrator: <input type="text" value="6"/>	OD: <input type="text"/>	Conc: <input type="text" value="**"/>		
Point 7:	Calibrator: <input type="text"/>	OD: <input type="text"/>	Conc: <input type="text"/>		
Point 8:	Calibrator: <input type="text"/>	OD: <input type="text"/>	Conc: <input type="text"/>		
Point 9:	Calibrator: <input type="text"/>	OD: <input type="text"/>	Conc: <input type="text"/>		
Point 10:	Calibrator: <input type="text"/>	OD: <input type="text"/>	Conc: <input type="text"/>		
<Point Cal. For		No. of Correction Points: <input type="text"/>	Use Master Curve: <input type="checkbox"/>	Lot Calibration: <input type="checkbox"/>	
Master Curve>					
Point-1:	Calibrator: <input type="text"/>	OD: <input type="text"/>	Conc: <input type="text"/>	Stability Reagent Blank: <input type="text" value="28"/> Day	<input type="text" value="0"/> Hour
Point-2:	Calibrator: <input type="text"/>	OD: <input type="text"/>	Conc: <input type="text"/>	Stability Calibration: <input type="text" value="28"/> Day	<input type="text" value="0"/> Hour
MB Type Factor: <input type="text"/>		1-Point Calibration Point: <input type="checkbox"/>	<input type="checkbox"/> with Conc-0		

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Systemreagens: B08179  
 Testnaam: CYS1G

Reagens ID: 228

General	LIH	ISE	Calculated Test	Range
<b>Test Name:</b> CYS1G <input type="button" value="Test No"/> <b>Type:</b> Serum <b>Operation:</b> Yes				
Sample Volume	<input type="text" value="2.0"/> µL	Dilution	<input type="text" value="0"/> µL	OD Limit
Pre-Dilution Rate	<input type="text" value="1"/>			Min. OD <input type="text"/> Max OD <input type="text"/>
Reagent Volume	R1 (R1-1) <input type="text" value="150"/> µL	Dilution	<input type="text" value="0"/> µL	Reagent OD Limit 1 <sup>st</sup> Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R1-2 <input type="text"/> µL	Dilution	<input type="text"/> µL	Last Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R2 (R2-1) <input type="text" value="30"/> µL	Dilution	<input type="text" value="10"/> µL	Analytical Measuring Range Low <input type="text" value="0.4"/> High <input type="text" value="7.8"/>
Common Reagent	Type <input type="text" value="None"/>	Name	<input type="text" value="None"/>	Correlation Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec	<input type="text" value="None"/> nm	Manufacturer Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Method	<input type="text" value="END"/>			
Reaction Slope	<input type="text" value="+"/> ▼			Onboard Stability Period <input type="text" value="60"/> Day <input type="text" value="0"/> Hour
Measuring Point-1	1st <input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check <input type="text" value="No"/> ▼
Measuring Point-2	1st <input type="text"/>	Last	<input type="text"/>	Lipemia <input type="text"/> ▼
Linearity Limit	<input type="text"/> %			Icterus <input type="text"/> ▼
Lag Time Check	<input type="text"/> ▼			Hemolysis <input type="text"/> ▼

General	LIH	ISE	Calculated Test	Range
<b>Test Name:</b> CYS1G <input type="button" value="Test No"/> <b>Type:</b> Serum				
Value/Flag	<input type="text" value="Value"/>	Level	Low <input type="text" value="-99999.99"/> High <input type="text" value="99999.99"/>	
Specific Ranges				
	From	To		
<input type="checkbox"/> 1:	Sex <input type="text" value="#"/> ▼ Year <input type="text" value="#"/> Month <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	Other Type <input type="text" value="None"/>	Low <input type="text" value="#"/> High <input type="text" value="#"/>
<input type="checkbox"/> 2:	<input type="text" value="#"/> ▼ <input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 3:	<input type="text" value="#"/> ▼ <input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 4:	<input type="text" value="#"/> ▼ <input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 5:	<input type="text" value="#"/> ▼ <input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 6:	<input type="text" value="#"/> ▼ <input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
7:	Standard demographics			<input type="text" value="#"/> <input type="text" value="#"/>
8:	Not within expected values			<input type="text" value="#"/> <input type="text" value="#"/>
Critical Limits	Low <input type="text" value="#"/>	High <input type="text" value="#"/>	Unit <input type="text" value="mg/L"/>	<input type="button" value="Select"/> Decimal Places <input type="text" value="2"/>

Calibrators	General	ISE
<b>Test Name:</b> CYS1G <input type="checkbox"/> Use Serum Cal. <b>Type:</b> Serum		
Calibration Type:	<input type="text" value="6AB"/> ▼	Formula: <input type="text" value="Spline"/> ▼
		Counts: <input type="text" value="2"/>
<Calibrator Parameters>		Slope Check <input type="text" value="+"/> ▼
		Allowable Range Check
		<input type="checkbox"/> Reagent Blank <input type="text"/>
		<input type="checkbox"/> Calibration <input type="text"/>
		Advanced Calibration
		Operation <input type="text" value="No"/> ▼
		Interval (RB) <input type="text"/> ▼
		Interval (ACAL) <input type="text"/> ▼
MB Type Factor	<input type="text"/>	1-Point Calibration Point <input type="text" value="None"/> ▼ <input type="checkbox"/> with Conc=0
		Stability
		Reagent Blank <input type="text" value="28"/> Day <input type="text" value="0"/> Hour
		Calibration <input type="text" value="28"/> Day <input type="text" value="0"/> Hour

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