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Systemreagenz: B08179 Reagenz-ID: 228

Specific Test Parameters										
General		LIH	ISE	Range						
Test Name:	CysC	<	>	Type:	Serum	Operation:	Yes			
Sample:	Volume	2	µL	Dilution	0	µL	Pre-Dilution Rate:	1		
Reagents:	R1 Volume	150	µL	Dilution	0	µL	Min OD		Max OD	
	R2 Volume	30	µL	Dilution	10	µL		L	H	
Wavelength:	Pri.	540	nm	Sec.		nm	Reagent OD limit:			
Method:	End Point						First L	-2.0	First H	2.0
Reaction slope:		+					Last L		Last H	
Measuring Point 1:	First	13		Last	27		Dynamic Range:	L	H	*
Measuring Point 2:	First			Last			Correlation Factor:	A	B	1.00 0.00
Linearity:			%				On-board stability period:	60		
No Lag Time:										

Specific Test Parameters									
General		LIH	ISE	Range					
Test Name:	CysC	<	>	Type:	Serum				
Value/Flag:	#	Level L:	#	Level H:	#				
Normal Ranges:	Age L		Age H		L		H		
	Sex	Year	Month	Year	Month				
o 1.	#	#	#	#	#	#	#	#	#
o 2.	#	#	#	#	#	#	#	#	#
o 3.	#	#	#	#	#	#	#	#	#
o 4.	#	#	#	#	#	#	#	#	#
o 5.	#	#	#	#	#	#	#	#	#
o 6.	#	#	#	#	#	#	#	#	#
7.	None Selected								
8.	Out of Range		L	H		#	#	#	#
Panic Value:	#	#	#	#	#	Unit:	mg/L	Decimal places:	#

Calibration Specific									
General		ISE							
Test Name:	CysC	<	>	Type:	Serum				
Calibration Type:	6AB	Formula:	Spline	Counts:	#	Process:	CONC		
Point 1:	Cal. No.	OD	CONC	Factor/OD-L	Factor/OD-H				
Point 2:	1		**						
Point 3:	2		**						
Point 4:	3		**						
Point 5:	4		**						
Point 6:	5		**						
Point 7:	6		**						
1-Point Cal. Point:		o	with CONC-0	Slope Check:	+	Advanced Calibration:	#		
MB Type Factor:			Calibration Stability Period:	28					

User defined

** Lot specific, see analytical value sheet included with calibrator kit

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Systemreagenz: B08179 Reagenz-ID: 228

Parameters		Specific Test Parameters		
General	LIH	ISE	Calculated Test	Range
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation: <input type="text" value="Yes"/>				
Sample Volume	<input type="text" value="2"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD <input type="text"/> Max.OD <input type="text"/>
Rgt. Volume	R1(R1-1) <input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit
				First Low <input type="text" value="-2.0"/> High <input type="text" value="2.0"/>
				Last Low <input type="text"/> High <input type="text"/>
	R2(R2-1) <input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Dynamic Range Low <input type="text"/> High <input type="text"/>
Common Rgt. Type	<input type="text"/>	Name	<input type="text"/>	Correlation Factor A <input type="text" value="1.00"/> B <input type="text" value="0.00"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec.	<input type="text"/>	Factor for Maker A <input type="text"/> B <input type="text"/>
Method	<input type="text" value="End Point"/>			
Reaction Slope	<input type="text" value="+"/> ∇			Onboard Stability Period <input type="text" value="60"/> Day <input type="text"/> Hour
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check <input type="text"/>
Measuring Point2 First	<input type="text"/>	Last	<input type="text"/>	Lipemia <input type="text"/>
Linearity Limit	<input type="text"/>			Icterus <input type="text"/>
Lag Time Check	<input type="text"/>			Hemolysis <input type="text"/>

Parameters		Specific Test Parameters		
General	LIH	ISE	Calculated Test	Range
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>				
Value/Flag:	<input type="text" value="#"/>	Low	<input type="text"/>	High <input type="text"/>
Level	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Specific Ranges:	From	To		Panic Value Low <input type="text"/> High <input type="text"/>
o 1.	Sex <input type="text"/>	Year <input type="text"/>	Month <input type="text"/>	Year <input type="text"/>
o 2.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
o 3.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
o 4.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
o 5.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
o 6.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
o 7.	No demographics			<input type="text"/>
o 8.	Not within expected values			<input type="text"/>
Unit	<input type="text" value="mg/L"/>	Decimal Places	<input type="text" value="#"/>	

Parameters		Calibration Parameters		
Calibrators	Calibration Specific	STAT Table Calibration		
General	ISE			
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> <input type="checkbox"/> Use Serum Cal.				
Calibration Type:	<input type="text" value="6AB"/>	Formula:	<input type="text" value="Spline"/>	Counts: <input type="text" value="#"/>
<Calibrator Parameters>				
Calibrator	OD	Conc	Factor Range Low High	Slope Check <input type="text" value="+"/>
Point 1:	<input type="text" value="1"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 2:	<input type="text" value="2"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 3:	<input type="text" value="3"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 4:	<input type="text" value="4"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 5:	<input type="text" value="5"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 6:	<input type="text" value="6"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 7:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 8:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 9:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 10:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<Point Cal. For	No. of Correction Points <input type="text"/>	Use Master Curve <input type="checkbox"/>	Lot Calibration <input type="checkbox"/>	
Master Curve>				
Calibrator	OD	Conc	Low High	Stability
Point-1	<input type="text"/>	<input type="text"/>	<input type="text"/>	Reagent Blank <input type="text" value="28"/> Day <input type="text"/> Hour
Point-2	<input type="text"/>	<input type="text"/>	<input type="text"/>	Calibration <input type="text" value="28"/> Day <input type="text"/> Hour
MB Type Factor:	<input type="text"/>	1-Point Calibration Point <input type="checkbox"/>	<input type="text"/>	<input type="checkbox"/> with Conc-0

User defined
 ** Lot specific, see analytical value sheet included with calibrator kit

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Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Cuvette: <input type="text" value=""/>					
<input type="checkbox"/> Use Serum Cal.					
Calibration Type: <input type="text" value="6AB"/> Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="#"/>					
<Calibrator Parameters>					
Point 1:	Calibrator	OD	Conc	Factor Range	Slope Check
Point 2:	<input type="text" value="1"/>	<input type="text" value=""/>	<input type="text" value="**"/>	Low High	<input type="text" value="+"/>
Point 3:	<input type="text" value="2"/>	<input type="text" value=""/>	<input type="text" value="**"/>		Allowance Range Check
Point 4:	<input type="text" value="3"/>	<input type="text" value=""/>	<input type="text" value="**"/>		<input type="checkbox"/> Reagent Blank
Point 5:	<input type="text" value="4"/>	<input type="text" value=""/>	<input type="text" value="**"/>		<input type="checkbox"/> Calibration
Point 6:	<input type="text" value="5"/>	<input type="text" value=""/>	<input type="text" value="**"/>		Advanced Calibration
Point 7:	<input type="text" value="6"/>	<input type="text" value=""/>	<input type="text" value="**"/>		Operation <input type="text" value="#"/>
Point 8:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>		Interval (RB/ACAL) <input type="text" value="#"/>
Point 9:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>		
Point 10:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>		
<Point Cal. For No. of Correction Points <input type="text" value=""/> Use Master Curve <input type="checkbox"/> Lot Calibration <input type="checkbox"/>					
Master Curve>					
Point-1	Calibrator	OD	Conc	Low High	Stability
Point-2	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>		Reagent Blank <input type="text" value="28"/> Day <input type="text" value="0"/> Hour
					Calibration <input type="text" value="28"/> Day <input type="text" value="0"/> Hour
MB Type Factor: <input type="text" value=""/> 1-Point Calibration Point <input type="checkbox"/> with Conc-0					

User defined

** Lot specific, see analytical value sheet included with calibrator kit