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Specific Test Parameters											
General		LIH	ISE	Range							
Test Name:		CysC		<	>	Type:	Serum		Operation:	Yes	
Sample:	Volume	2	µL	Dilution	0	µL	Pre-Dilution Rate:	1			
Reagents:	R1 Volume	150	µL	Dilution	0	µL	Min OD	Max OD			
	R2 Volume	30	µL	Dilution	10	µL	L	H			
Wavelength:		Pri.	540	Sec.			Reagent OD limit:				
Method:		End Point				First L		-2.0	First H		2.0
Reaction slope:		+				Last L		Last H			
Measuring Point 1:		First	13	Last	27		Dynamic Range:				
Measuring Point 2:		First			Last		L		H		*
Linearity:						Correlation Factor:		A	1.00	B	0.00
No Lag Time:						On-board stability period:		60			

Specific Test Parameters										
General		LIH	ISE	Range						
Test Name:		CysC		<	>	Type:	Serum			
Value/Flag:		#	Level L:	#	Level H:	#				
Normal Ranges:		Age L		Age H						
	Sex	Year	Month	Year	Month	L	H			
o 1.	#	#	#	#	#	#	#			
o 2.	#	#	#	#	#	#	#			
o 3.	#	#	#	#	#	#	#			
o 4.	#	#	#	#	#	#	#			
o 5.	#	#	#	#	#	#	#			
o 6.	#	#	#	#	#	#	#			
7. None Selected						#	#			
8. Out of Range		L		H		#	#			
Panic Value:		#	L	#	H	Unit:	mg/L		Decimal places: #	

Calibration Specific											
General		ISE									
Test Name:		CysC		<	>	Type:	Serum				
Calibration Type:		6AB		Formula:	Spline		Counts:	#		Process:	CONC
Point 1:	Cal. No.	1	OD	CONC	Factor/OD-L	Factor/OD-H					
Point 2:	2			**							
Point 3:	3			**							
Point 4:	4			**							
Point 5:	5			**							
Point 6:	6			**							
Point 7:				**							
1-Point Cal. Point:		#	with CONC-0	Slope Check:	+		Advanced Calibration:	#			
MB Type Factor:				Calibration Stability Period:		28					

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Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation <input type="text" value="Yes"/>					
Sample Volume	<input type="text" value="2"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD	Max.OD
Rgt. Volume	R1(R1-1) <input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit	
				First	Low
				Last	Low
				High	<input type="text" value="2.0"/>
				High	<input type="text" value="2.0"/>
R2(R2-1)	<input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Dynamic Range Low	
				High	<input type="text" value="0.00"/>
Common Rgt. Type	<input type="text" value=""/>	Name	<input type="text" value=""/>	Correlation Factor A	<input type="text" value="1.00"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec.	<input type="text" value=""/>	Factor for Maker A	<input type="text" value=""/>
Method	<input type="text" value="End Point"/>				
Reaction Slope	<input type="text" value="+"/> ▾	Onboard Stability Period		<input type="text" value="60"/> Day	<input type="text" value=""/>
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check <input type="text" value=""/>	
Measuring Point2 First	<input type="text" value=""/>	Last	<input type="text" value=""/>	Lipemia <input type="text" value=""/>	
Linearity Limit	<input type="text" value=""/>	Icterus <input type="text" value=""/>			
Lag Time Check	<input type="text" value=""/>	Hemolysis <input type="text" value=""/>			

Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>					
Value/Flag:	<input type="text" value="#"/>	Low	<input type="text" value="#"/>	High <input type="text" value="#"/>	
Level	<input type="text" value="#"/>	From	<input type="text" value="#"/>	To	<input type="text" value="#"/>
Specific Ranges:					
	Sex	Year	Month	Year	Month
o 1.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
o 2.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
o 3.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
o 4.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
o 5.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
o 6.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
o 7.	No demographics				
o 8.	Not within expected values				
Unit	<input type="text" value="mg/L"/>	Decimal Places	<input type="text" value="#"/>	Panic Value	
				Low	High
				<input type="text" value="#"/>	<input type="text" value="#"/>

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> <input type="checkbox"/> Use Serum Cal.					
Calibration Type: <input type="text" value="6AB"/> Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="#"/>					
<Calibrator Parameters>					
	Calibrator	OD	Conc	Low	High
Point 1:	<input type="text" value="1"/>	<input type="text" value=""/>	<input type="text" value="**"/>	<input type="text" value=""/>	<input type="text" value=""/>
Point 2:	<input type="text" value="2"/>	<input type="text" value=""/>	<input type="text" value="**"/>	<input type="text" value=""/>	<input type="text" value=""/>
Point 3:	<input type="text" value="3"/>	<input type="text" value=""/>	<input type="text" value="**"/>	<input type="text" value=""/>	<input type="text" value=""/>
Point 4:	<input type="text" value="4"/>	<input type="text" value=""/>	<input type="text" value="**"/>	<input type="text" value=""/>	<input type="text" value=""/>
Point 5:	<input type="text" value="5"/>	<input type="text" value=""/>	<input type="text" value="**"/>	<input type="text" value=""/>	<input type="text" value=""/>
Point 6:	<input type="text" value="6"/>	<input type="text" value=""/>	<input type="text" value="**"/>	<input type="text" value=""/>	<input type="text" value=""/>
Point 7:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
Point 8:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
Point 9:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
Point 10:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
Slope Check <input type="text" value="+"/>					
Allowance Range Check					
o Reagent Blank <input type="text" value=""/>					
o Calibration <input type="text" value=""/>					
Advanced Calibration					
Operation <input type="text" value="#"/>					
Interval (RB/ACAL) <input type="text" value="#"/>					
<Point Cal. For No. of Correction Points <input type="text" value=""/>					
Use Master Curve <input type="checkbox"/> Lot Calibration <input type="checkbox"/>					
Master Curve>					
	Calibrator	OD	Conc	Low	High
Point-1	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
Point-2	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
Stability					
Reagent Blank <input type="text" value="28"/> Day <input type="text" value=""/>					
Calibration <input type="text" value="28"/> Day <input type="text" value=""/>					
MB Type Factor: <input type="text" value=""/>					
1-Point Calibration Point <input type="checkbox"/> with Conc-0 <input type="checkbox"/>					

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