Cistatinas C

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Sistemos reagentas: B08179 Reagento ID: 228

Specific Test Parameters									
General		LIH	ISE	Range					
Test Name:	CysC	<	>	Type:	Serum	Operation:	Yes		
Sample:	Volume	2	μL	Dilution	0	μL	Pre-Dilution Rate:	1	
Reagents:	R1 Volume	150	μL	Dilution	0	μL	Min OD		Max OD
	R2 Volume	30	μL	Dilution	10	μL	L		H
Wavelength:	Pri.	540	nm	Sec.		nm	Reagent OD limit:		
Method:	End Point			First L	-2.0	First H	2.0	Last L	Last H
Reaction slope:		+		Dynamic Range:	L	*	H	*	
Measuring Point 1:	First	13		Last	27		Correlation Factor:	A	1.00
Measuring Point 2:	First			Last			B	0.00	
Linearity:			%	On-board stability period:	60				
No Lag Time:									

Specific Test Parameters									
General		LIH	ISE	Range					
Test Name:	CysC	<	>	Type:	Serum				
Value/Flag:	#	Level L:	#	Level H:	#				
Normal Ranges:	Age L		Age H		L		H		
	Sex	Year	Month	Year	Month				
o 1.	#	#	#	#	#	#	#	#	#
o 2.	#	#	#	#	#	#	#	#	#
o 3.	#	#	#	#	#	#	#	#	#
o 4.	#	#	#	#	#	#	#	#	#
o 5.	#	#	#	#	#	#	#	#	#
o 6.	#	#	#	#	#	#	#	#	#
7.	None Selected								
8.	Out of Range		L	H					
Panic Value:	#	#	#	Unit:	mg/L	Decimal places:	#		

Calibration Specific									
General		ISE							
Test Name:	CysC	<	>	Type:	Serum				
Calibration Type:	6AB	Formula:	Spline	Counts:	#	Process:	CONC		
Point 1:	Cal. No.	OD	CONC	Factor/OD-L	Factor/OD-H				
Point 2:	1		**						
Point 3:	2		**						
Point 4:	3		**						
Point 5:	4		**						
Point 6:	5		**						
Point 7:	6		**						
1-Point Cal. Point:		o	with CONC-0	Slope Check:	+	Advanced Calibration:	#		
MB Type Factor:		Calibration Stability Period:	28						

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Sistemos reagentas: B08179 Reagento ID: 228

Parameters		Specific Test Parameters		
General	LIH	ISE	Calculated Test	Range
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation <input type="text" value="Yes"/>				
Sample Volume	<input type="text" value="2"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD <input type="text"/> Max.OD <input type="text"/>
Rgt. Volume	R1(R1-1) <input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit
				First Low <input type="text" value="-2.0"/> High <input type="text" value="2.0"/>
				Last Low <input type="text"/> High <input type="text"/>
	R2(R2-1) <input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Dynamic Range Low <input type="text"/> High <input type="text"/>
Common Rgt. Type	<input type="text"/>	Name	<input type="text"/>	Correlation Factor A <input type="text" value="1.00"/> B <input type="text" value="0.00"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec.	<input type="text"/>	Factor for Maker A <input type="text"/> B <input type="text"/>
Method	<input type="text" value="End Point"/>			
Reaction Slope	<input type="text" value="+"/> ∇			Onboard Stability Period <input type="text" value="60"/> Day <input type="text"/> Hour
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check <input type="text"/>
Measuring Point2 First	<input type="text"/>	Last	<input type="text"/>	Lipemia <input type="text"/>
Linearity Limit	<input type="text"/>			Icterus <input type="text"/>
Lag Time Check	<input type="text"/>			Hemolysis <input type="text"/>

Parameters		Specific Test Parameters		
General	LIH	ISE	Calculated Test	Range
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>				
Value/Flag:	<input type="text" value="#"/>	Low	<input type="text"/>	High <input type="text"/>
Level	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Specific Ranges:	From	To		Panic Value Low <input type="text"/> High <input type="text"/>
1. Sex	<input type="text"/>	<input type="text"/>	<input type="text"/>	Low <input type="text"/> High <input type="text"/>
2. Year	<input type="text"/>	<input type="text"/>	<input type="text"/>	Low <input type="text"/> High <input type="text"/>
3. Month	<input type="text"/>	<input type="text"/>	<input type="text"/>	Low <input type="text"/> High <input type="text"/>
4. Year	<input type="text"/>	<input type="text"/>	<input type="text"/>	Low <input type="text"/> High <input type="text"/>
5. Month	<input type="text"/>	<input type="text"/>	<input type="text"/>	Low <input type="text"/> High <input type="text"/>
6. Year	<input type="text"/>	<input type="text"/>	<input type="text"/>	Low <input type="text"/> High <input type="text"/>
7. No demographics	<input type="text"/>	<input type="text"/>	<input type="text"/>	Low <input type="text"/> High <input type="text"/>
8. Not within expected values	<input type="text"/>	<input type="text"/>	<input type="text"/>	Low <input type="text"/> High <input type="text"/>
Unit	<input type="text" value="mg/L"/>	Decimal Places	<input type="text" value="#"/>	

Parameters		Calibration Parameters		
Calibrators	Calibration Specific	STAT Table Calibration		
General	ISE			
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> <input type="checkbox"/> Use Serum Cal.				
Calibration Type:	<input type="text" value="6AB"/>	Formula:	<input type="text" value="Spline"/>	Counts: <input type="text" value="#"/>
<Calibrator Parameters>				
Calibrator	OD	Conc	Factor Range Low High	Slope Check <input type="text" value="+"/> ∇
Point 1:	<input type="text" value="1"/>	<input type="text"/>	<input type="text"/>	Allowance Range Check <input type="checkbox"/> Reagent Blank <input type="text"/> <input type="checkbox"/> Calibration <input type="text"/> Advanced Calibration Operation <input type="text" value="#"/> ∇ Interval (RB/ACAL) <input type="text" value="#"/> ∇
Point 2:	<input type="text" value="2"/>	<input type="text"/>	<input type="text"/>	
Point 3:	<input type="text" value="3"/>	<input type="text"/>	<input type="text"/>	
Point 4:	<input type="text" value="4"/>	<input type="text"/>	<input type="text"/>	
Point 5:	<input type="text" value="5"/>	<input type="text"/>	<input type="text"/>	
Point 6:	<input type="text" value="6"/>	<input type="text"/>	<input type="text"/>	
Point 7:	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 8:	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 9:	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 10:	<input type="text"/>	<input type="text"/>	<input type="text"/>	
<Point Cal. For	No. of Correction Points <input type="text"/>	Use Master Curve <input type="checkbox"/>	Lot Calibration <input type="checkbox"/>	
Master Curve>				
Calibrator	OD	Conc	Low High	Stability
Point-1	<input type="text"/>	<input type="text"/>	<input type="text"/>	Reagent Blank <input type="text" value="28"/> Day <input type="text"/> Hour
Point-2	<input type="text"/>	<input type="text"/>	<input type="text"/>	Calibration <input type="text" value="28"/> Day <input type="text"/> Hour
MB Type Factor:	<input type="text"/>	1-Point Calibration Point <input type="checkbox"/>	<input type="text"/>	<input type="checkbox"/> with Conc-0

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Sistemos reagentas: B08179 Reagento ID: 228

Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation <input type="text" value="Yes"/>					
Sample Volume	<input type="text" value="2"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>	Diluent Bottle	<input type="text" value="Outside"/>	Min.OD	<input type="text" value=""/>
Rgt. Volume	<input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit	
R1-2	<input type="text" value=""/>	Dilution	<input type="text" value=""/>	First	Low <input type="text" value="-2.0"/> High <input type="text" value="2.0"/>
R2(R2-1)	<input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Last	Low <input type="text" value=""/>
Common Rgt. Type	<input type="text" value=""/>	Name	<input type="text" value=""/>	Dynamic Range Low	<input type="text" value="*"/> High <input type="text" value="*"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec.	<input type="text" value=""/>	Correlation Factor A	<input type="text" value="1.00"/> B <input type="text" value="0.00"/>
Method	<input type="text" value="End Point"/>	Factor for Maker A			
Reaction Slope	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	Onboard Stability Period	<input type="text" value="60"/> Day <input type="text" value=""/>
Measuring Point1 First		LIH Influence Check <input type="text" value=""/>			
Measuring Point2 First		Lipemia <input type="text" value=""/>			
Linearity Limit	<input type="text" value=""/>	Icterus <input type="text" value=""/>			
Lag Time Check	<input type="text" value=""/>	Hemolysis <input type="text" value=""/>			

Parameters		Specific Test Parameters					
General	LIH	ISE	Calculated Test	Range			
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>							
Value/Flag: <input type="text" value="#"/>							
Level Low <input type="text" value="#"/> High <input type="text" value="#"/>							
Specific Ranges: From To Low High							
1.	Sex	Year	Month	Year	Month	Low	High
2.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
3.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
4.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
5.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
6.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
7.	Standard demographics					<input type="text" value="#"/>	<input type="text" value="#"/>
8.	Not within expected values					<input type="text" value="#"/>	<input type="text" value="#"/>
Panic Value Low <input type="text" value=""/> High <input type="text" value=""/> Unit <input type="text" value="mg/L"/> Decimal Places <input type="text" value="#"/>							

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Cuvette: <input type="text" value=""/>					
<input type="checkbox"/> Use Serum Cal.					
Calibration Type: <input type="text" value="6AB"/> Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="#"/>					
<Calibrator Parameters>					
Point 1:	Calibrator	OD	Conc	Factor	Range
Point 2:	<input type="text" value="1"/>	<input type="text" value=""/>	<input type="text" value="**"/>	Low	High
Point 3:	<input type="text" value="2"/>	<input type="text" value=""/>	<input type="text" value="**"/>		
Point 4:	<input type="text" value="3"/>	<input type="text" value=""/>	<input type="text" value="**"/>		
Point 5:	<input type="text" value="4"/>	<input type="text" value=""/>	<input type="text" value="**"/>		
Point 6:	<input type="text" value="5"/>	<input type="text" value=""/>	<input type="text" value="**"/>		
Point 7:	<input type="text" value="6"/>	<input type="text" value=""/>	<input type="text" value="**"/>		
Point 8:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>		
Point 9:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>		
Point 10:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>		
Slope Check <input type="text" value="+"/> Allowance Range Check					
<input type="checkbox"/> Reagent Blank <input type="text" value=""/>					
<input type="checkbox"/> Calibration <input type="text" value=""/>					
Advanced Calibration Operation <input type="text" value="#"/>					
Interval (RB/ACAL) <input type="text" value="#"/>					
<Point Cal. For No. of Correction Points <input type="text" value=""/> Use Master Curve <input type="checkbox"/> Lot Calibration <input type="checkbox"/>					
Master Curve>					
Point-1	Calibrator	OD	Conc	Low	High
Point-2	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>		
Stability Reagent Blank <input type="text" value="28"/> Day <input type="text" value="0"/> Hour					
Calibration <input type="text" value="28"/> Day <input type="text" value="0"/> Hour					
MB Type Factor: <input type="text" value=""/> 1-Point Calibration Point <input type="checkbox"/> with Conc-0					

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