Cistatina C

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Specific Test Parameters											
General		LIH	ISE	Range							
Test Name:		CysC		<	>	Type:	Serum		Operation:	Yes	
Sample:	Volume	2	µL	Dilution	0	µL	Pre-Dilution Rate:	1			
Reagents:	R1 Volume	150	µL	Dilution	0	µL	Min OD	Max OD			
	R2 Volume	30	µL	Dilution	10	µL	L	H			
Wavelength:	Pri.	540	∇	Sec.		∇	Reagent OD limit:				
Method:	End Point		∇	First L	-2.0		First H	2.0			
Reaction slope:	+		∇	Last L			Last H				
Measuring Point 1:	First	13		Last	27		Dynamic Range:				
Measuring Point 2:	First			Last			L	0.4		H	7.8
Linearity:			%	Correlation Factor:			A	1.00		B	0.00
No Lag Time:			∇	On-board stability period:	60						

Specific Test Parameters										
General		LIH	ISE	Range						
Test Name:		CysC		<	>	Type:	Serum			
Value/Flag:	#	∇	Level L:	#		Level H:	#			
Normal Ranges:	Sex	Age L	Year	Month	Age H	Year	Month	L	H	
<input type="checkbox"/>	1.	#	∇	#	#	#	#	#	#	#
<input type="checkbox"/>	2.	#	∇	#	#	#	#	#	#	#
<input type="checkbox"/>	3.	#	∇	#	#	#	#	#	#	#
<input type="checkbox"/>	4.	#	∇	#	#	#	#	#	#	#
<input type="checkbox"/>	5.	#	∇	#	#	#	#	#	#	#
<input type="checkbox"/>	6.	#	∇	#	#	#	#	#	#	#
	7.	None Selected								
	8.	Out of Range		L	H		#	#	#	#
Panic Value:	#		#	Unit:	mg/L		Decimal places:	#		

Calibration Specific											
General		ISE									
Test Name:		CysC		<	>	Type:	Serum				
Calibration Type:	6AB		∇	Formula:	Spline		∇	Counts:	#		
				Process:	CONC						
Point 1:	Cal. No.	1	OD	CONC	Factor/OD-L	Factor/OD-H					
Point 2:	2			**							
Point 3:	3			**							
Point 4:	4			**							
Point 5:	5			**							
Point 6:	6			**							
Point 7:				**							
1-Point Cal. Point:	<input type="checkbox"/>	with CONC-0	Slope Check:	+		∇	Advanced Calibration:	#			
MB Type Factor:			Calibration Stability Period:	28							

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Parameters		Specific Test Parameters		
General	LIH	ISE	Calculated Test	Range
Test Name: <input type="text" value="CysC"/> <input type="button" value="<"/> <input type="button" value=">"/> Type: <input type="text" value="Serum"/> Operation: <input type="text" value="Yes"/>				
Sample Volume	<input type="text" value="2"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD <input type="text"/> Max.OD <input type="text"/>
Rgt. Volume	R1(R1-1) <input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit
				First Low <input type="text" value="-2.0"/> High <input type="text" value="2.0"/>
				Last Low <input type="text"/> High <input type="text"/>
R2(R2-1)	<input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	
Common Rgt. Type	<input type="text"/>	Name	<input type="text"/>	Dynamic Range Low <input type="text" value="0.4"/> High <input type="text" value="7.8"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec. <input type="text"/>	<input type="text"/>	Correlation Factor A <input type="text" value="1.00"/> B <input type="text" value="0.00"/>
Method	<input type="text" value="End Point"/>			Factor for Maker A <input type="text"/>
Reaction Slope	<input type="text" value="+"/> ∇			Onboard Stability Period <input type="text" value="60"/> Day <input type="text"/> Hour
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check <input type="text"/>
Measuring Point2 First	<input type="text"/>	Last	<input type="text"/>	Lipemia <input type="text"/>
Linearity Limit	<input type="text"/>			Icterus <input type="text"/>
Lag Time Check	<input type="text"/>			Hemolysis <input type="text"/>

Parameters		Specific Test Parameters		
General	LIH	ISE	Calculated Test	Range
Test Name: <input type="text" value="CysC"/> <input type="button" value="<"/> <input type="button" value=">"/> Type: <input type="text" value="Serum"/>				
Value/Flag:	<input type="text" value="#"/>	Low	<input type="text"/>	High <input type="text"/>
Level	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Specific Ranges:	From	To		
<input type="checkbox"/> 1.	Sex <input type="text"/>	Year <input type="text"/>	Month <input type="text"/>	Year <input type="text"/>
<input type="checkbox"/> 2.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 3.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 4.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 5.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 6.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 7.	No demographics			<input type="text"/>
<input type="checkbox"/> 8.	Not within expected values			<input type="text"/>
Unit	<input type="text" value="mg/L"/>	Decimal Places	<input type="text" value="#"/>	

Parameters		Calibration Parameters		
Calibrators	Calibration Specific	STAT Table Calibration		
General	ISE			
Test Name: <input type="text" value="CysC"/> <input type="button" value="<"/> <input type="button" value=">"/> Type: <input type="text" value="Serum"/> <input type="checkbox"/> Use Serum Cal.				
Calibration Type:	<input type="text" value="6AB"/>	Formula:	<input type="text" value="Spline"/>	Counts: <input type="text" value="#"/>
<Calibrator Parameters>				
Calibrator	OD	Conc	Factor Range	Slope Check <input type="text" value="+"/> ∇
Point 1:	<input type="text" value="1"/>	<input type="text"/>	<input type="text"/>	Allowance Range Check
Point 2:	<input type="text" value="2"/>	<input type="text"/>	<input type="text"/>	
Point 3:	<input type="text" value="3"/>	<input type="text"/>	<input type="text"/>	
Point 4:	<input type="text" value="4"/>	<input type="text"/>	<input type="text"/>	
Point 5:	<input type="text" value="5"/>	<input type="text"/>	<input type="text"/>	
Point 6:	<input type="text" value="6"/>	<input type="text"/>	<input type="text"/>	
Point 7:	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 8:	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 9:	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 10:	<input type="text"/>	<input type="text"/>	<input type="text"/>	
<input type="checkbox"/> Reagent Blank <input type="text"/>				
<input type="checkbox"/> Calibration <input type="text"/>				
Advanced Calibration Operation <input type="text" value="#"/> ∇				
Interval (RB/ACAL) <input type="text" value="#"/> ∇				
<Point Cal. For Master Curve>				
No. of Correction Points	<input type="text"/>	Use Master Curve	<input type="checkbox"/>	<input type="checkbox"/> Lot Calibration
Calibrator	OD	Conc	Low	High
Point-1	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point-2	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Stability	Reagent Blank	<input type="text" value="28"/>	Day	<input type="text"/>
	Calibration	<input type="text" value="28"/>	Day	<input type="text"/>
MB Type Factor:	<input type="text"/>	1-Point Calibration Point	<input type="checkbox"/>	<input type="checkbox"/> with Conc-0

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Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation: <input type="text" value="Yes"/>					
Sample Volume	<input type="text" value="2"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>	Diluent Bottle	<input type="text" value="Outside"/>	Min.OD	<input type="text"/>
Rgt. Volume	<input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit	
R1(R1-1)	<input type="text"/>	Dilution	<input type="text"/>	First	Low <input type="text" value="-2.0"/> High <input type="text" value="2.0"/>
R1-2	<input type="text"/>	Dilution	<input type="text"/>	Last	Low <input type="text"/>
R2(R2-1)	<input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Dynamic Range Low	<input type="text" value="0.4"/> High <input type="text" value="7.8"/>
Common Rgt. Type	<input type="text"/>	Name	<input type="text"/>	Correlation Factor A	<input type="text" value="1.00"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec.	<input type="text"/>	Factor for Maker A	<input type="text"/>
Method	<input type="text" value="End Point"/>			Factor for Maker B	<input type="text" value="0.00"/>
Reaction Slope	<input type="text"/>			Onboard Stability Period	<input type="text" value="60"/> Day <input type="text"/> Hour
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check	<input type="text"/>
Measuring Point2 First	<input type="text"/>	Last	<input type="text"/>	Lipemia	<input type="text"/>
Linearity Limit	<input type="text"/>			Icterus	<input type="text"/>
Lag Time Check	<input type="text"/>			Hemolysis	<input type="text"/>

Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>					
Value/Flag:	<input type="text" value="#"/>				
Level		Level	Low <input type="text" value="#"/>	High <input type="text" value="#"/>	
Specific Ranges:	From	To	Low	High	
<input type="checkbox"/> 1.	Sex <input type="text" value="#"/>	Year <input type="text" value="#"/>	Month <input type="text" value="#"/>	Year <input type="text" value="#"/>	Month <input type="text" value="#"/>
<input type="checkbox"/> 2.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 3.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 4.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 5.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 6.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 7.	Standard demographics				
<input type="checkbox"/> 8.	Not within expected values				
Panic Value	Low <input type="text"/>	High <input type="text"/>	Unit <input type="text" value="mg/L"/>	Decimal Places <input type="text" value="#"/>	

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Cuvette: <input type="text"/>					
<input type="checkbox"/> Use Serum Cal.					
Calibration Type: <input type="text" value="6AB"/> Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="#"/>					
<Calibrator Parameters>					
Calibrator	OD	Conc	Factor Range		Slope Check <input type="text" value="+"/>
			Low	High	
Point 1:	<input type="text" value="1"/>	<input type="text" value="**"/>			Allowance Range Check <input type="checkbox"/> Reagent Blank <input type="text"/> <input type="checkbox"/> Calibration <input type="text"/> Advanced Calibration Operation <input type="text" value="#"/> Interval (RB/ACAL) <input type="text" value="#"/>
Point 2:	<input type="text" value="2"/>	<input type="text" value="**"/>			
Point 3:	<input type="text" value="3"/>	<input type="text" value="**"/>			
Point 4:	<input type="text" value="4"/>	<input type="text" value="**"/>			
Point 5:	<input type="text" value="5"/>	<input type="text" value="**"/>			
Point 6:	<input type="text" value="6"/>	<input type="text" value="**"/>			
Point 7:	<input type="text"/>				
Point 8:	<input type="text"/>				
Point 9:	<input type="text"/>				
Point 10:	<input type="text"/>				
<Point Cal. For No. of Correction Points <input type="text"/> Use Master Curve <input type="text"/> Lot Calibration <input type="checkbox"/>					
Master Curve>					
Calibrator	OD	Conc	OD Range		Stability
			Low	High	Reagent Blank <input type="text" value="28"/> Day <input type="text" value="0"/> Hour
Point-1	<input type="text"/>	<input type="text"/>			Calibration <input type="text" value="28"/> Day <input type="text" value="0"/> Hour
Point-2	<input type="text"/>	<input type="text"/>			
MB Type Factor: <input type="text"/> 1-Point Calibration Point <input type="text"/> with Conc-0 <input type="checkbox"/>					

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General	LIH	ISE	Calculated Test	Range
Test Name: CYS1G <input type="button" value="Test No"/> Type: Serum <input type="button" value="Operation"/> Yes <input type="button" value="▼"/>				
Sample Volume	<input type="text" value="2.0"/> µL	Dilution	<input type="text" value="0"/> µL	OD Limit
Pre-Dilution Rate	<input type="text" value="1"/> ▼			Min. OD <input type="text"/> Max OD <input type="text"/>
Reagent Volume	R1 (R1-1) <input type="text" value="150"/> µL	Dilution	<input type="text" value="0"/> µL	Reagent OD Limit 1 st Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R1-2 <input type="text"/> µL	Dilution	<input type="text"/> µL	Last Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R2 (R2-1) <input type="text" value="30"/> µL	Dilution	<input type="text" value="10"/> µL	Analytical Measuring Range Low <input type="text" value="0.4"/> High <input type="text" value="7.8"/>
Common Reagent	Type <input type="text" value="None"/>	Name	<input type="text" value="None"/>	Correlation Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec	<input type="text" value="None"/> nm	Manufacturer Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Method	<input type="text" value="END"/> ▼			
Reaction Slope	<input type="text" value="+"/> ▼			Onboard Stability Period <input type="text" value="60"/> Day <input type="text" value="0"/> Hour
Measuring Point-1	1st <input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check <input type="text" value="No"/> ▼
Measuring Point-2	1st <input type="text"/>	Last	<input type="text"/>	Lipemia <input type="text"/> ▼
Linearity Limit	<input type="text"/> %			Icterus <input type="text"/> ▼
Lag Time Check	<input type="text"/> ▼			Hemolysis <input type="text"/> ▼

General	LIH	ISE	Calculated Test	Range
Test Name: CYS1G <input type="button" value="Test No"/> Type: Serum <input type="button" value="▼"/>				
Value/Flag	<input type="text" value="Value"/>	Level	Low <input type="text" value="-99999.99"/>	High <input type="text" value="99999.99"/>
Specific Ranges				
	From	To		
	Sex	Year	Month	Other Type
<input type="checkbox"/> 1:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
<input type="checkbox"/> 2:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
<input type="checkbox"/> 3:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
<input type="checkbox"/> 4:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
<input type="checkbox"/> 5:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
<input type="checkbox"/> 6:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
7:	Standard demographics			<input type="text" value="#"/>
8:	Not within expected values			<input type="text" value="#"/>
Critical Limits	Low <input type="text" value="#"/>	High <input type="text" value="#"/>	Unit <input type="text" value="mg/L"/>	<input type="button" value="Select"/> Decimal Places <input type="text" value="2"/>

Calibrators	General	ISE
Test Name: CYS1G <input type="checkbox" value="Use Serum Cal."/>		
Calibration Type:	<input type="text" value="6AB"/> ▼	Formula: <input type="text" value="Spline"/> ▼
		Counts: <input type="text" value="2"/>
<Calibrator Parameters>		Slope Check <input type="text" value="+"/> ▼
		Allowable Range Check
		<input type="checkbox"/> Reagent Blank <input type="text"/>
		<input type="checkbox"/> Calibration <input type="text"/>
		Advanced Calibration
		Operation <input type="text" value="No"/> ▼
		Interval (RB) <input type="text"/> ▼
		Interval (ACAL) <input type="text"/> ▼
MB Type Factor	<input type="text"/>	1-Point Calibration Point <input type="text" value="None"/> ▼
		<input type="checkbox"/> with Conc=0
		Stability
		Reagent Blank <input type="text" value="28"/> Day <input type="text" value="0"/> Hour
		Calibration <input type="text" value="28"/> Day <input type="text" value="0"/> Hour

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