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Specific Test Parameters										
General		LIH	ISE	Range						
Test Name:	CysC	<	>	Type:	Serum	Operation:	Yes			
Sample:	Volume	2	μL	Dilution	0	μL	Pre-Dilution Rate:	1		
Reagents:	R1 Volume	150	μL	Dilution	0	μL	Min OD	Max OD		
	R2 Volume	30	μL	Dilution	10	μL	L	H		
Wavelength:	Pri.	540	▽	Sec.		▽	Reagent OD limit:			
Method:	End Point		▽	First L	-2.0	First H	2.0			
Reaction slope:		+	▽	Last L		Last H				
Measuring Point 1:	First	13		Last	27		Dynamic Range:	L 0.4 H 7.8		
Measuring Point 2:	First			Last			Correlation Factor:	A 1.00 B 0.00		
Linearity:			%	On-board stability period:						60
No Lag Time:			▽							

Specific Test Parameters										
General		LIH	ISE	Range						
Test Name:	CysC	<	>	Type:	Serum					
Value/Flag:	#	▽	Level L:	#	Level H:	#				
Normal Ranges:	Sex	Age L	Year	Month	Age H	Year	Month	L	H	
<input type="checkbox"/>	1.	#	▽	#	#	#	#	#	#	
<input type="checkbox"/>	2.	#	▽	#	#	#	#	#	#	
<input type="checkbox"/>	3.	#	▽	#	#	#	#	#	#	
<input type="checkbox"/>	4.	#	▽	#	#	#	#	#	#	
<input type="checkbox"/>	5.	#	▽	#	#	#	#	#	#	
<input type="checkbox"/>	6.	#	▽	#	#	#	#	#	#	
	7.	None Selected								
	8.	Out of Range								
Panic Value:		#	L	#	H	#	Unit:	mg/L	Decimal places:	#

Calibration Specific										
General		ISE								
Test Name:	CysC	<	>	Type:	Serum					
Calibration Type:	6AB	▽	Formula:	Spline	▽	Counts:	#	Process:	CONC	▽
Point 1:	Cal. No.	1	OD	CONC	Factor/OD-L	Factor/OD-H				
Point 2:	2			**						
Point 3:	3			**						
Point 4:	4			**						
Point 5:	5			**						
Point 6:	6			**						
Point 7:				**						
1-Point Cal. Point:	<input type="checkbox"/>	with CONC-0	Slope Check:	+	▽	Advanced Calibration:	#	▽		
MB Type Factor:			Calibration Stability Period:	28						

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Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation: <input type="text" value="Yes"/>					
Sample Volume	<input type="text" value="2"/> μ L	Dilution	<input type="text" value="0"/> μ L	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD	<input type="text"/>
Rgt. Volume	R1(R1-1) <input type="text" value="150"/> μ L	Dilution	<input type="text" value="0"/> μ L	Reagent OD Limit	
				First	Low <input type="text" value="-2.0"/> High <input type="text" value="2.0"/>
				Last	Low <input type="text"/> High <input type="text"/>
	R2(R2-1) <input type="text" value="30"/> μ L	Dilution	<input type="text" value="10"/> μ L	Dynamic Range Low <input type="text" value="0.4"/> High <input type="text" value="7.8"/>	
Common Rgt. Type	<input type="text"/>	Name	<input type="text"/>	Correlation Factor A	<input type="text" value="1.00"/>
Wavelength	Pri <input type="text" value="540"/> ∇ nm	Sec.	<input type="text"/>	Factor for Maker A	B <input type="text" value="0.00"/>
Method	<input type="text" value="End Point"/>				
Reaction Slope	<input type="text" value="+"/>			Onboard Stability Period	<input type="text" value="60"/> Day <input type="text"/> Hour
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check	<input type="text"/>
Measuring Point2 First	<input type="text"/>	Last	<input type="text"/>	Lipemia	<input type="text"/>
Linearity Limit	<input type="text"/>			Icterus	<input type="text"/>
Lag Time Check	<input type="text"/>			Hemolysis	<input type="text"/>

Parameters		Specific Test Parameters					
General	LIH	ISE	Calculated Test	Range			
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>							
Value/Flag:	<input type="text" value="#"/>	Low	<input type="text"/>	High	<input type="text"/>		
Level	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		
Specific Ranges:	From	To	Low	High	Panic Value		
					Low <input type="text"/> High <input type="text"/>		
<input type="checkbox"/> 1.	Sex <input type="text"/>	Year <input type="text"/>	Month <input type="text"/>	Year <input type="text"/>	Month <input type="text"/>	Low <input type="text"/>	High <input type="text"/>
<input type="checkbox"/> 2.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 3.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 4.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 5.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 6.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 7.	No demographics		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 8.	Not within expected values		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Unit	<input type="text" value="mg/L"/>		Decimal Places	<input type="text" value="#"/>			

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> <input type="checkbox"/> Use Serum Cal.					
Calibration Type: <input type="text" value="6AB"/> Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="#"/>					
<Calibrator Parameters>					
Calibrator	OD	Conc	Low	High	Slope Check <input type="text" value="+"/>
Point 1:	<input type="text" value="1"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Allowance Range Check <input type="checkbox"/> Reagent Blank <input type="text"/> <input type="checkbox"/> Calibration <input type="text"/> Advanced Calibration Operation <input type="text" value="#"/> Interval (RB/ACAL) <input type="text" value="#"/>
Point 2:	<input type="text" value="2"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 3:	<input type="text" value="3"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 4:	<input type="text" value="4"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 5:	<input type="text" value="5"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 6:	<input type="text" value="6"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 7:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 8:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 9:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 10:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
<Point Cal. For No. of Correction Points <input type="text"/> Use Master Curve <input type="checkbox"/> Lot Calibration <input type="checkbox"/>					
<Master Curve>					
Calibrator	OD	Conc	Low	High	Stability
Point-1	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Reagent Blank <input type="text" value="28"/> Day <input type="text"/>
Point-2	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Calibration <input type="text" value="28"/> Day <input type="text"/>
MB Type Factor: <input type="text"/> 1-Point Calibration Point <input type="checkbox"/> with Conc-0 <input type="checkbox"/>					

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General	LIH	ISE	Calculated Test	Range
Test Name: CYS1G Test No Type: Serum Operation: Yes				
Sample Volume	<input type="text" value="2.0"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit
Pre-Dilution Rate	<input type="text" value="1"/>			Min. OD <input type="text"/> Max OD <input type="text"/>
Reagent Volume	R1 (R1-1) <input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit 1 st Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R1-2 <input type="text"/> μL	Dilution	<input type="text"/> μL	Last Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R2 (R2-1) <input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Analytical Measuring Range Low <input type="text" value="0.4"/> High <input type="text" value="7.8"/>
Common Reagent	Type <input type="text" value="None"/>	Name	<input type="text" value="None"/>	Correlation Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
	Wavelength Pri <input type="text" value="540"/> nm	Sec <input type="text" value="None"/> nm	Manufacturer Factor A <input type="text" value="1"/> B <input type="text" value="0"/>	
Method	<input type="text" value="END"/>			
Reaction Slope	<input type="text" value="+"/>		Onboard Stability Period <input type="text" value="60"/> Day <input type="text" value="0"/> Hour	
Measuring Point-1	1st <input type="text" value="13"/>	Last <input type="text" value="27"/>	LIH Influence Check <input type="text" value="No"/>	Lipemia <input type="text"/>
Measuring Point-2	1st <input type="text"/>	Last <input type="text"/>	Icterus <input type="text"/>	Hemolysis <input type="text"/>
Linearity Limit	<input type="text"/> %			
Lag Time Check	<input type="text"/>			

General	LIH	ISE	Calculated Test	Range
Test Name: CYS1G Test No Type: Serum				
Value/Flag	<input type="text" value="Value"/>	Level	Low <input type="text" value="-99999.99"/> High <input type="text" value="99999.99"/>	
Specific Ranges				
	From	To	Other Type	Low High
<input type="checkbox"/> 1:	Sex <input type="text" value="#"/> Year <input type="text" value="#"/> Month <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 2:	<input type="text" value="#"/> <input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 3:	<input type="text" value="#"/> <input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 4:	<input type="text" value="#"/> <input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 5:	<input type="text" value="#"/> <input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 6:	<input type="text" value="#"/> <input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
7:	Standard demographics			<input type="text" value="#"/> <input type="text" value="#"/>
8:	Not within expected values			<input type="text" value="#"/> <input type="text" value="#"/>
Critical Limits	Low <input type="text" value="#"/>	High <input type="text" value="#"/>	Unit <input type="text" value="mg/L"/>	Select <input type="text"/> Decimal Places <input type="text" value="2"/>

Calibrators	General	ISE		
Test Name: CYS1G Type: Serum				
<input type="checkbox"/> Use Serum Cal.				
Calibration Type:	<input type="text" value="6AB"/>	Formula:	<input type="text" value="Spline"/>	Counts: <input type="text" value="2"/>
<Calibrator Parameters>				Slope Check <input type="text" value="+"/>
	Calibrator	OD	Conc	Range Low High
Point-1	CYSC Calibrator Level 1		†	-2.0000 2.0000
Point-2	CYSC Calibrator Level 2		†	-2.0000 2.0000
Point-3	CYSC Calibrator Level 3		†	-2.0000 2.0000
Point-4	CYSC Calibrator Level 4		†	-2.0000 2.0000
Point-5	CYSC Calibrator Level 5		†	-2.0000 2.0000
Point-6	CYSC Calibrator Level 6		†	-2.0000 2.0000
Point-7				
MB Type Factor	<input type="text"/>	1-Point Calibration Point	<input type="text" value="None"/>	<input type="checkbox"/> with Conc-0
				Stability Reagent Blank <input type="text" value="28"/> Day <input type="text" value="0"/> Hour
				Calibration <input type="text" value="28"/> Day <input type="text" value="0"/> Hour

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