Cystatin C

Folheto informativo para Imunoensaio Cistatina C nos Sistemas AU da Beckman Coulter® (AU5800, AU680, AU480, AU2700, DxC 700 AU)

REF B08179

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Reagente de Sistema: B08179 ID de Reagente: 228

Specific Test Parameters										
General		LIH	ISE	Range						
Test Name:	CysC	<	>	Type:	Serum	Operation:	Yes			
Sample:	Volume	2	µL	Dilution	0	µL	Pre-Dilution Rate:	1		
Reagents:	R1 Volume	150	µL	Dilution	0	µL	Min OD	Max OD		
	R2 Volume	30	µL	Dilution	10	µL	L	H		
Wavelength:	Pri.	540	∇	Sec.		∇	Reagent OD limit:			
Method:	End Point		∇	First L	-2.0	First H	2.0			
Reaction slope:	+		∇	Last L		Last H				
Measuring Point 1:	First	13		Last	27		Dynamic Range:			
Measuring Point 2:	First			Last			L	0.4	H	7.8
Linearity:			%			Correlation Factor:				
No Lag Time:			∇			A	1.00	B	0.00	
				On-board stability period:		60				

Specific Test Parameters										
General		LIH	ISE	Range						
Test Name:	CysC	<	>	Type:	Serum					
Value/Flag:	#	∇	Level L:	#	Level H:	#				
Normal Ranges:	Age L	Age H			L	H				
	Sex	Year	Month	Year	Month	L	H			
<input type="checkbox"/>	1.	#	#	#	#	#	#			
<input type="checkbox"/>	2.	#	#	#	#	#	#			
<input type="checkbox"/>	3.	#	#	#	#	#	#			
<input type="checkbox"/>	4.	#	#	#	#	#	#			
<input type="checkbox"/>	5.	#	#	#	#	#	#			
<input type="checkbox"/>	6.	#	#	#	#	#	#			
7. None Selected										
8. Out of Range										
Panic Value:	L		#	H		#	Unit:	mg/L	Decimal places:	#

Calibration Specific										
General		ISE								
Test Name:	CysC	<	>	Type:	Serum					
Calibration Type:	6AB	∇	Formula:	Spline	∇	Counts:	#	Process:	CONC	∇
Point 1:	Cal. No.	1	OD	CONC	**	Factor/OD-L	Factor/OD-H			
Point 2:	2			**						
Point 3:	3			**						
Point 4:	4			**						
Point 5:	5			**						
Point 6:	6			**						
Point 7:				**						
1-Point Cal. Point:	<input type="checkbox"/>	with CONC-0	Slope Check:	+	∇	Advanced Calibration:	#			
MB Type Factor:			Calibration Stability Period:	28						

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Reagente de Sistema: B08179 ID de Reagente: 228

Parameters		Specific Test Parameters		
General	LIH	ISE	Calculated Test	Range
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation: <input type="text" value="Yes"/>				
Sample Volume	<input type="text" value="2"/> μ L	Dilution	<input type="text" value="0"/> μ L	OD Limit
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD <input type="text"/> Max.OD <input type="text"/>
Rgt. Volume	R1(R1-1) <input type="text" value="150"/> μ L	Dilution	<input type="text" value="0"/> μ L	Reagent OD Limit
				First Low <input type="text" value="-2.0"/> High <input type="text" value="2.0"/>
				Last Low <input type="text"/> High <input type="text"/>
R2(R2-1)	<input type="text" value="30"/> μ L	Dilution	<input type="text" value="10"/> μ L	
Common Rgt. Type	<input type="text"/>	Name	<input type="text"/>	Dynamic Range Low <input type="text" value="0.4"/> High <input type="text" value="7.8"/>
Wavelength	Pri <input type="text" value="540"/> ν nm	Sec.	<input type="text"/>	Correlation Factor A <input type="text" value="1.00"/> B <input type="text" value="0.00"/>
Method	<input type="text" value="End Point"/>			Factor for Maker A <input type="text"/>
Reaction Slope	<input type="text" value="+"/> ν			Onboard Stability Period <input type="text" value="60"/> Day <input type="text"/> Hour
Measuring Point1	First <input type="text" value="13"/> Last <input type="text"/>	Last	<input type="text" value="27"/>	LIH Influence Check <input type="text"/>
Measuring Point2	First <input type="text"/>	Last	<input type="text"/>	Lipemia <input type="text"/>
Linearity Limit	<input type="text"/>			Icterus <input type="text"/>
Lag Time Check	<input type="text"/>			Hemolysis <input type="text"/>

Parameters		Specific Test Parameters		
General	LIH	ISE	Calculated Test	Range
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>				
Value/Flag:	<input type="text" value="#"/>	Low	<input type="text"/>	High <input type="text"/>
Level				
Specific Ranges:	From	To		Panic Value
				Low <input type="text"/> High <input type="text"/>
<input type="checkbox"/> 1.	Sex <input type="text"/>	Year <input type="text"/>	Month <input type="text"/>	Year <input type="text"/>
<input type="checkbox"/> 2.				
<input type="checkbox"/> 3.				
<input type="checkbox"/> 4.				
<input type="checkbox"/> 5.				
<input type="checkbox"/> 6.				
<input type="checkbox"/> 7.	No demographics			
<input type="checkbox"/> 8.	Not within expected values			
Unit	<input type="text" value="mg/L"/>	Decimal Places	<input type="text" value="#"/>	

Parameters		Calibration Parameters		
Calibrators	Calibration Specific	STAT Table Calibration		
General	ISE			
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> <input type="checkbox"/> Use Serum Cal.				
Calibration Type:	<input type="text" value="6AB"/>	Formula:	<input type="text" value="Spline"/>	Counts: <input type="text" value="#"/>
<Calibrator Parameters>				
Calibrator	OD	Conc	Factor Range	Slope Check <input type="text" value="+"/>
Point 1:	<input type="text" value="1"/>	<input type="text"/>	<input type="text"/>	
Point 2:	<input type="text" value="2"/>	<input type="text"/>	<input type="text"/>	Allowance Range Check
Point 3:	<input type="text" value="3"/>	<input type="text"/>	<input type="text"/>	
Point 4:	<input type="text" value="4"/>	<input type="text"/>	<input type="text"/>	<input type="checkbox"/> Reagent Blank <input type="text"/>
Point 5:	<input type="text" value="5"/>	<input type="text"/>	<input type="text"/>	<input type="checkbox"/> Calibration <input type="text"/>
Point 6:	<input type="text" value="6"/>	<input type="text"/>	<input type="text"/>	Advanced Calibration
Point 7:	<input type="text"/>	<input type="text"/>	<input type="text"/>	Operation <input type="text" value="#"/>
Point 8:	<input type="text"/>	<input type="text"/>	<input type="text"/>	Interval (RB/ACAL) <input type="text" value="#"/>
Point 9:	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 10:	<input type="text"/>	<input type="text"/>	<input type="text"/>	
<Point Cal. For Master Curve>	No. of Correction Points <input type="text"/>	Use Master Curve <input type="checkbox"/>	Lot Calibration <input type="checkbox"/>	
Calibrator	OD	Conc	Factor Range	Stability
Point-1	<input type="text"/>	<input type="text"/>	<input type="text"/>	Reagent Blank <input type="text" value="28"/> Day <input type="text"/>
Point-2	<input type="text"/>	<input type="text"/>	<input type="text"/>	Calibration <input type="text" value="28"/> Day <input type="text"/>
MB Type Factor:	<input type="text"/>	1-Point Calibration Point <input type="checkbox"/>	with Conc-0 <input type="checkbox"/>	

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ID de Reagente: 228

General	LIH	ISE	Calculated Test	Range
Test Name: CYS1G <input type="button" value="Test No"/> Type: Serum <input type="button" value="Operation"/> Yes				
Sample Volume	<input type="text" value="2.0"/> µL	Dilution	<input type="text" value="0"/> µL	OD Limit
Pre-Dilution Rate	<input type="text" value="1"/>			Min. OD <input type="text"/> Max OD <input type="text"/>
Reagent Volume	R1 (R1-1) <input type="text" value="150"/> µL	Dilution	<input type="text" value="0"/> µL	Reagent OD Limit 1 st Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R1-2 <input type="text"/> µL	Dilution	<input type="text"/> µL	Last Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R2 (R2-1) <input type="text" value="30"/> µL	Dilution	<input type="text" value="10"/> µL	Analytical Measuring Range Low <input type="text" value="0.4"/> High <input type="text" value="7.8"/>
Common Reagent	Type <input type="text" value="None"/>	Name	<input type="text" value="None"/>	Correlation Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec	<input type="text" value="None"/> nm	Manufacturer Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Method	<input type="text" value="END"/>			
Reaction Slope	<input type="text" value="+"/> ▾			Onboard Stability Period <input type="text" value="60"/> Day <input type="text" value="0"/> Hour
Measuring Point-1	1st <input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check <input type="text" value="No"/> ▾
Measuring Point-2	1st <input type="text"/>	Last	<input type="text"/>	Lipemia <input type="text"/> ▾
Linearity Limit	<input type="text"/> %			Icterus <input type="text"/> ▾
Lag Time Check	<input type="text"/> ▾			Hemolysis <input type="text"/> ▾

General	LIH	ISE	Calculated Test	Range
Test Name: CYS1G <input type="button" value="Test No"/> Type: Serum <input type="button" value="Operation"/>				
Value/Flag	<input type="text" value="Value"/>	Level	Low <input type="text" value="-99999.99"/> High <input type="text" value="99999.99"/>	
Specific Ranges				
	Sex	Year	Month	Year
	From	To		
<input type="checkbox"/> 1:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 2:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 3:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 4:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 5:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 6:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
7:	Standard demographics			<input type="text" value="#"/>
8:	Not within expected values			<input type="text" value="#"/>
Critical Limits	Low <input type="text" value="#"/>	High <input type="text" value="#"/>	Unit <input type="text" value="mg/L"/>	<input type="button" value="Select"/> Decimal Places <input type="text" value="2"/>

Calibrators	General	ISE		
Test Name: CYS1G <input type="button" value="Use Serum Cal."/>				
Calibration Type:	<input type="text" value="6AB"/>	Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="2"/>		
<Calibrator Parameters>				
	Calibrator	OD		
		Conc		
		Range		
		Low		
		High		
Point-1	CYSC Calibrator Level 1 ▾	†	-2.0000	2.0000
Point-2	CYSC Calibrator Level 2 ▾	†	-2.0000	2.0000
Point-3	CYSC Calibrator Level 3 ▾	†	-2.0000	2.0000
Point-4	CYSC Calibrator Level 4 ▾	†	-2.0000	2.0000
Point-5	CYSC Calibrator Level 5 ▾	†	-2.0000	2.0000
Point-6	CYSC Calibrator Level 6 ▾	†	-2.0000	2.0000
Point-7				
MB Type Factor	<input type="text"/>	1-Point Calibration Point <input type="text" value="None"/> ▾	<input type="checkbox"/> with Conc=0	Slope Check <input type="text" value="+"/> ▾
				Allowable Range Check
				<input type="checkbox"/> Reagent Blank <input type="text"/>
				<input type="checkbox"/> Calibration <input type="text"/>
				Advanced Calibration
				Operation <input type="text" value="No"/> ▾
				Interval (RB) <input type="text"/> ▾
				Interval (ACAL) <input type="text"/> ▾
				Stability
				Reagent Blank <input type="text" value="28"/> Day <input type="text" value="0"/> Hour
				Calibration <input type="text" value="28"/> Day <input type="text" value="0"/> Hour

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