Bipacksedel för Cystatin C-immunanlys i Beckman Coulter® AU-system (AU5800, AU680, AU480, AU2700, DxC 700 AU)

REF B08179

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Systemreagens: B08179 Reagens-ID: 228

Specific Test Parameters											
General		LIH	ISE	Range							
Test Name:	CysC	<	>	Type:	Serum	Operation:	Yes				
Sample:	Volume	2	µL	Dilution	0	µL	Pre-Dilution Rate:	1			
Reagents:	R1 Volume	150	µL	Dilution	0	µL	Min OD	Max OD			
	R2 Volume	30	µL	Dilution	10	µL	L	H			
Wavelength:	Pri.	540	∇	Sec.		∇	Reagent OD limit:				
Method:	End Point						First L	-2.0	First H	2.0	
Reaction slope:	+						Last L		Last H		
Measuring Point 1:	First	13		Last	27		Dynamic Range:	L	0.4	H	7.8
Measuring Point 2:	First			Last			Correlation Factor:				
Linearity:							A	1.00	B	0.00	
No Lag Time:							On-board stability period:	60			

Specific Test Parameters										
General		LIH	ISE	Range						
Test Name:	CysC	<	>	Type:	Serum					
Value/Flag:	#	∇	Level L:	#	Level H:	#				
Normal Ranges:	Age L	Age H			L		H			
	Sex	Year	Month	Year	Month					
<input type="checkbox"/>	1.	#	#	#	#	#	#	#	#	
<input type="checkbox"/>	2.	#	#	#	#	#	#	#	#	
<input type="checkbox"/>	3.	#	#	#	#	#	#	#	#	
<input type="checkbox"/>	4.	#	#	#	#	#	#	#	#	
<input type="checkbox"/>	5.	#	#	#	#	#	#	#	#	
<input type="checkbox"/>	6.	#	#	#	#	#	#	#	#	
7. None Selected										
8. Out of Range										
Panic Value:	L		#	H		#	Unit:	mg/L	Decimal places:	#

Calibration Specific									
General		ISE							
Test Name:	CysC	<	>	Type:	Serum				
Calibration Type:	6AB	Formula:	Spline	Counts:	#	Process:	CONC		
Point 1:	Cal. No.	OD	CONC	Factor/OD-L	Factor/OD-H				
Point 2:	1		**						
Point 3:	2		**						
Point 4:	3		**						
Point 5:	4		**						
Point 6:	5		**						
Point 7:	6		**						
Point 7:	7		**						
1-Point Cal. Point:	<input type="checkbox"/>	with CONC-0	Slope Check:	+	Advanced Calibration:	#			
MB Type Factor:		Calibration Stability Period:	28						

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Systemreagens: B08179 Reagens-ID: 228

Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation: <input type="text" value="Yes"/>					
Sample Volume	<input type="text" value="2"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD	<input type="text"/>
Rgt. Volume	R1(R1-1) <input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit	
				First	Low <input type="text" value="-2.0"/> High <input type="text" value="2.0"/>
				Last	Low <input type="text"/> High <input type="text"/>
R2(R2-1)	<input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Dynamic Range Low	<input type="text" value="0.4"/> High <input type="text" value="7.8"/>
Common Rgt. Type	<input type="text"/>	Name	<input type="text"/>	Correlation Factor A	<input type="text" value="1.00"/> B <input type="text" value="0.00"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec.	<input type="text"/>	Factor for Maker A	<input type="text"/>
Method	<input type="text" value="End Point"/>				
Reaction Slope	<input type="text" value="+"/> ∇			Onboard Stability Period	<input type="text" value="60"/> Day <input type="text"/> Hour
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check	<input type="text"/>
Measuring Point2 First	<input type="text"/>	Last	<input type="text"/>	Lipemia	<input type="text"/>
Linearity Limit	<input type="text"/>			Icterus	<input type="text"/>
Lag Time Check	<input type="text"/>			Hemolysis	<input type="text"/>

Parameters		Specific Test Parameters					
General	LIH	ISE	Calculated Test	Range			
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>							
Value/Flag:	<input type="text" value="#"/>	Low	<input type="text"/>	High	<input type="text"/>		
Level		<input type="text" value="#"/>		<input type="text" value="#"/>			
Specific Ranges:	From	To			Panic Value		
	Sex	Year	Month	Year	Month	Low	High
<input type="checkbox"/> 1.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 2.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 3.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 4.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 5.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 6.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
7.	No demographics					<input type="text" value="#"/>	<input type="text" value="#"/>
8.	Not within expected values					<input type="text" value="#"/>	<input type="text" value="#"/>
Unit	<input type="text" value="mg/L"/>		Decimal Places	<input type="text" value="#"/>			

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> <input type="checkbox"/> Use Serum Cal.					
Calibration Type: <input type="text" value="6AB"/> Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="#"/>					
<Calibrator Parameters>					
Calibrator	OD	Conc	Low	High	Slope Check
Point 1:	<input type="text" value="1"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="+"/>
Point 2:	<input type="text" value="2"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Allowance Range Check
Point 3:	<input type="text" value="3"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="checkbox"/> Reagent Blank
Point 4:	<input type="text" value="4"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="checkbox"/> Calibration
Point 5:	<input type="text" value="5"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Advanced Calibration
Point 6:	<input type="text" value="6"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Operation
Point 7:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="#"/>
Point 8:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Interval (RB/ACAL)
Point 9:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="#"/>
Point 10:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
<Point Cal. For Master Curve>					
No. of Correction Points	<input type="text"/>	Use Master Curve	<input type="checkbox"/>	Lot Calibration	<input type="checkbox"/>
Calibrator	OD	Conc	Low	High	Stability
Point-1	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Reagent Blank
Point-2	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Calibration
MB Type Factor:	<input type="text"/>	1-Point Calibration Point	<input type="checkbox"/>	with Conc-0	<input type="checkbox"/>

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Systemreagens: B08179
Testnamn: CYS1G

Reagens-ID: 228

General	LIH	ISE	Calculated Test	Range
Test Name: CYS1G <input type="button" value="Test No"/> Type: Serum <input type="button" value="Operation"/> Yes				
Sample Volume	<input type="text" value="2.0"/> µL	Dilution	<input type="text" value="0"/> µL	OD Limit
Pre-Dilution Rate	<input type="text" value="1"/>			Min. OD <input type="text"/> Max OD <input type="text"/>
Reagent Volume	R1 (R1-1) <input type="text" value="150"/> µL	Dilution	<input type="text" value="0"/> µL	Reagent OD Limit 1 st Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R1-2 <input type="text"/> µL	Dilution	<input type="text"/> µL	Last Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R2 (R2-1) <input type="text" value="30"/> µL	Dilution	<input type="text" value="10"/> µL	Analytical Measuring Range Low <input type="text" value="0.4"/> High <input type="text" value="7.8"/>
Common Reagent	Type <input type="text" value="None"/>	Name	<input type="text" value="None"/>	Correlation Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec	<input type="text" value="None"/> nm	Manufacturer Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Method	<input type="text" value="END"/>			
Reaction Slope	<input type="text" value="+"/> ▾			Onboard Stability Period <input type="text" value="60"/> Day <input type="text" value="0"/> Hour
Measuring Point-1	1st <input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check <input type="text" value="No"/> ▾
Measuring Point-2	1st <input type="text"/>	Last	<input type="text"/>	Lipemia <input type="text"/> ▾
Linearity Limit	<input type="text"/> %			Icterus <input type="text"/> ▾
Lag Time Check	<input type="text"/> ▾			Hemolysis <input type="text"/> ▾

General	LIH	ISE	Calculated Test	Range
Test Name: CYS1G <input type="button" value="Test No"/> Type: Serum <input type="button" value="Operation"/>				
Value/Flag	<input type="text" value="Value"/>	Level	Low <input type="text" value="-99999.99"/> High <input type="text" value="99999.99"/>	
Specific Ranges				
	From	To		
	Sex	Year	Month	Other Type
<input type="checkbox"/> 1:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
<input type="checkbox"/> 2:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
<input type="checkbox"/> 3:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
<input type="checkbox"/> 4:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
<input type="checkbox"/> 5:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
<input type="checkbox"/> 6:	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="None"/>
7:	Standard demographics			<input type="text" value="#"/>
8:	Not within expected values			<input type="text" value="#"/>
Critical Limits	Low <input type="text" value="#"/>	High <input type="text" value="#"/>	Unit <input type="text" value="mg/L"/>	<input type="button" value="Select"/> Decimal Places <input type="text" value="2"/>

Calibrators	General	ISE		
Test Name: CYS1G <input type="button" value="Use Serum Cal."/>				
Calibration Type:	<input type="text" value="6AB"/>	Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="2"/>		
<Calibrator Parameters>		Slope Check <input type="text" value="+"/> ▾		
	Calibrator	OD		
		Conc		
		Range		
		Low		
		High		
Point-1	CYSC Calibrator Level 1 ▾	†	-2.0000	2.0000
Point-2	CYSC Calibrator Level 2 ▾	†	-2.0000	2.0000
Point-3	CYSC Calibrator Level 3 ▾	†	-2.0000	2.0000
Point-4	CYSC Calibrator Level 4 ▾	†	-2.0000	2.0000
Point-5	CYSC Calibrator Level 5 ▾	†	-2.0000	2.0000
Point-6	CYSC Calibrator Level 6 ▾	†	-2.0000	2.0000
Point-7				
MB Type Factor	<input type="text"/>	1-Point Calibration Point <input type="text" value="None"/> ▾	<input type="checkbox"/> with Conc=0	Interval (RB) <input type="text"/> ▾
				Interval (ACAL) <input type="text"/> ▾
				Stability
				Reagent Blank <input type="text" value="28"/> Day <input type="text" value="0"/> Hour
				Calibration <input type="text" value="28"/> Day <input type="text" value="0"/> Hour

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